"Cardiopulmonary and metabolic stress during ICAROS Cloud training"

conducted at the Faculty of Sports Science by Anna Dmitrieva, Ruhr- University Bochum, Bochum, Germany (2021)

<u>Aim</u>

To compare differences between exercising on the ICAROS Cloud and exercising on a mat in terms of cardiopulmonary and metabolic parameters.

Findings

- Training on the ICAROS Cloud leads to an increase in cardiopulmonary and metabolic parameters compared to training on a training mat.
- The energy expenditure on the ICAROS Cloud is about 5% higher than on training on a mat.
- The lactate concentration after training on the ICAROS Cloud is 7 % higher than training on the mat. Similarly, the heart rate is also about 1 % higher on The ICAROS Cloud.
- Exercising on the ICAROS Cloud is more fun than exercising on a training mat.
- The findings clearly show that despite the tendency to be more demanding, there is a higher willingness to train on the ICAROS Cloud than the mat.







ICAROS

Methods

- Cross Sectional examination of exercising on the cloud vs exercising on a training mat.
- ICAROS Cloud and Training Mat are the independent variables
- All participants performed both the trainings ordered randomly and the below mentioned variables were measured.

Subjects

Participants	Mean Age	Mean Height	Mean Weight	
N = 16 (5 females, 11 males)	24 ± 1.5 years	176 ± 6.8 cm	71.3 ± 10.3 kg	

Independent Variables	Dependent Variables (Measured Variables)
 ICAROS Cloud – Unstable training surface Training Mat – Stable training surface 	 Energy Expenditure Heart Rate (HR) Lactate Concentration (MetaMax 3B Spirometer) Rate of Perceived Extraction (RPE) PACES Questionnaire





Study Design

Tasks	Duration	Content	Data Collection
1. Warm-up	1. 5 minutes	1. Bicycle Ergometer	1. Resting Lactate & HR
2. Exercise - Mat	2. 12 minutes	2. Exercise Execution	2. HR, Spirometry, RPE
3. Break	3. 30 minutes	3. Sitting	3. Lactate & HR (0, 1, 3 & 5 min)
4. Warm-up	4. 5 minutes	4. Bicycle Ergometer	4. Resting Lactate & HR
5. Exercise - Cloud	5. 12 minutes	5. Exercise Execution	5. HR, Spirometry, RPE
6. End of Load	6. 5 minutes	6. Sitting	6. Lactate & HR (0, 1, 3 & 5 min)



RESULTS

• Energy Expenditure (kcal/17min)

Energy expenditure within the exercise time and five-minute post load was about 5% higher for exercising on the ICAROS Cloud compared to exercising on the training mat albeit the effect size being small (Wilcoxon Test (p. = 0.005/d.= 0.9/95% CI [1.5; 6.4])).

- Mean values :
 - ICAROS Cloud : 83 ± 15 kcal /17 min
 - Training Mat : 79 ±16 kcal / 17 min
- Heart Rate (bpm)

T-tests for dependent samples also reveal no significance (p. = 0.209/95% CI [-1.7; 7.1]).

- Mean values :
 - ICAROS Cloud : 156 ± 16 bpm
 - Training Mat : 154 ± 18 bpm
- Lactate Concentration (mmol/l)

Lactate concentration in the blood is more while exercising on the ICAROS Cloud (T-test (p. = 0.314/95% Cl [-0.4; 1.2])).

- Mean values :
 - ICAROS Cloud : 5.8 ± 1.9 mmol/l
 - Training Mat : 5.4 ± 2.1 mmol/l







RESULTS

- <u>Rate of Perceived Exertion (Borg Scale)</u>
 - The t-test for dependent samples does not show significance (p. = 0.909/95% CI [- 0.7; 0.6])
 - Mean values :
 - ICAROS Cloud : 12.9 ± 1.7
 - Training Mat : 12.9 ± 2.0

• PACES Questionnaire

- Training on the ICAROS CLOUD is significantly better than training on the training mat in terms of fun factor(T-test (p. = 0.000/95% CI [0.3; 0.8]), with Cohen's d = 1.2.
- Mean Values :
 - ICAROS Cloud : 3.5 ± 0.6
 - Training Mat : 2.9 ± 0.4

Mean ± standard		Paired differences				
	deviation	Mean value	Mean value Standard deviation of the mean		95% confidence interval of the difference	
			ule mean	At	Upper	Sig. (2- sideways
RPE_CLOUD	12,9 ± 1,7	-,03571	1,22835	-,69026	,61883	,909
RPE_MATTE	12,9 ± 2,0					



ICAROS

Conclusions

- Training on the ICAROS Cloud leads to an increase in cardiopulmonary and metabolic parameters compared to training on a training mat as is evident in the absolute values of all the measured variables even though the differences were not significant.
- The energy metabolism on the ICAROS Cloud is about 5% higher than on training on a mat. Since the trainings on both the ICAROS Cloud and Mat mirror each other, the higher calorie consumption can directly be attributed to the effectiveness of ICAROS Cloud.
- The lactate concentration after training on the ICAROS Cloud is 7 % higher than training on the mat. Similarly, the heart rate is also about 1 % higher on The ICAROS Cloud.
- In terms of 'fun factor', there is a clear and significant difference between training on the Cloud and on the mat as it is more fun to train on the ICAROS Cloud.
- The findings show that despite the tendency to be more demanding, there is a higher willingness to train on the ICAROS Cloud than the mat. This would also imply a higher willingness to train on the ICAROS Cloud and therefore a potential to increase the intensity as well as the quantity of training sessions.